

















## NONABSORBABLE SURGICAL SUTURES

### DURACARE\* Monofilament Polypropylene

Virtues	:	Monofilament suture with excellent tensile strength and reasonable elasticity. Its pliable smooth surface makes it excellent for easy passage.
Material	:	Synthetic, Polypropylene Suture. Blue coloured for easy identification during surgery.
Tissue Reaction	:	Minimal tissue reaction. Does not support infection. No tissue rupture.
Applications	:	General closure, Cardiovascular, Ophthalmic, Orthopaedic, Plastic and Micro Surgeries.
Choice	:	Loop, Needled and pre-cut sterile sutures conforming to USP standards.
Caution	:	Avoid damage to the surface of suture.



NEEDLE LENGTH	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE IN USP				
			6-0(0.7)	5-0(1)	4-0(1.5)	3-0(2)	2-0(3)
 12 mm	3/8 Circle Reverse Cutting	45 cm	FS 878				
		70 cm		FS 880			
 16 mm	3/8 Circle Round Bodied	70 cm		FS 881			
 16 mm	3/8 Circle Cutting	70 cm			FS 870		
 16 mm	1/2 Circle Round Bodied	70 cm			FS 849		
 17 mm, 17 mm Double Needle	1/2 Circle Taper cutting	70 cm			FS 835		
 20 mm	1/2 Circle Round Bodied	70 cm				FS 887	
 22 mm	3/8 Circle Cutting	80 cm				FS 018	FS 017
 25 mm	1/2 Circle Round Bodied	90 cm				FS 825	
 25 mm, 25 mm Double Needle	1/2 Circle Round Bodied	70 cm				FS 822	

NEEDLE LENGTH	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE IN USP		
			2-0(3)	1-0(3.5)	1(4)
	1/2 Circle Taper cutting	70 cm		FS 805	
 Double Needle	1/2 Circle Taper Cutting	90 cm	FS 877		
 Double Needle	1/2 Circle Round Bodied	90 cm		FS 824	
	1/2 Circle Round Bodied	90 cm	FS 841	FS 842	
	1/2 Circle Round Bodied	90 cm		FS 830	FS 840
	1/2 Circle Round Bodied	100 cm			FS 843
	1/2 Circle Reverse Cutting	100 cm			FS 883