
















NONABSORBABLE SURGICAL SUTURES

| LINEX* | Monofilament Polyamide |
|-----------------|---|
| Virtues | : Being monofilament, Linex is remarkably smooth and gives excellent knot security. Possesses good tensile strength. Easy removal with no tissue adherence. |
| Material | : Polyamide 6 suture (Nylon). Blue/Black coloured for easy identification in the operating field. |
| Tissue Reaction | : Linex Suture is inert like stainless steel. Maintains tensile strength indefinitely in tissues. Does not support infection. |
| Applications | : General closure, Skin and Plastic Surgery. Also used for Ophthalmic and Micro Surgery. |
| Choice | : Loop, Needled and pre-cut sterile sutures conforming to USP standards. |
| Caution | : Not recommended for attachment of artificial prostheses in Cardiovascular surgery. |



| NEEDLE LENGTH | NEEDLE DESCRIPTION | SUTURE LENGTH | SUTURE SIZE IN USP | | | | | | |
|---|-----------------------------------|---------------|--------------------|----------|---------|----------|----------|----------|--|
| | | | 10-0 | 6-0(0.7) | 5-0(1) | 4-0(1.5) | 3-0(2) | 2-0(3) | |
|  6 mm 6 mm Double needle | 3/8 Circle spatulated micro point | 38 cm | FS 3719 | | | | | | |
|  | 3/8 Circle Reverse Cutting | 70 cm | | FS 3320 | | | | | |
|  | 3/8 Circle Reverse Cutting | 70 cm | | | FS 3317 | | | | |
|  | 3/8 Circle Cutting | 35 cm | | | FS 3315 | | | | |
| | | 70 cm | | | | FS 3318 | FS 3321 | | |
|  | 3/8 Circle Reverse Cutting | 100 cm | | | | | FS 3328 | FS 3329 | |
| | | 45 cm | | | | | FS 3328S | FS 3329S | |

| NEEDLE LENGTH | NEEDLE DESCRIPTION | SUTURE LENGTH | SUTURE SIZE IN USP | | |
|---|----------------------------|---------------|--------------------|----------|---------|
| | | | 2-0(3) | 1-0(3.5) | 1(4) |
|  | 3/8 Circle Cutting | 100 cm | FS 3236 | FS 3239 | |
|  | 1/2 Circle Reverse cutting | 100 cm | | FS 3337 | |
| | | 45 cm | | FS 3337S | |
|  | 1/2 Circle Round Bodied | 150 cm (Loop) | | FS 3340 | |
|  | 1/2 Circle Reverse cutting | 100 cm | | FS 3339 | |
|  | 1/2 Circle Round Bodied | 100 cm | | FS 3346 | FS 3347 |

| NEEDLE LENGTH | NEEDLE DESCRIPTION | SUTURE LENGTH | SUTURE SIZE IN USP | | |
|---|----------------------------|---------------|--------------------|----------|---------|
| | | | 2-0(3) | 1-0(3.5) | 1(4) |
|  45 mm | 3/8 Circle Reverse Cutting | 45 cm | FS 3336S | | |
| | | 100 cm | FS 3336 | FS 3335 | FS 3334 |
|  45 mm Heavy | 1/2 Circle Reverse Cutting | 100 cm | | | FS 3338 |
|  50 mm Heavy | 1/2 Circle Round Bodied | 150 cm (Loop) | | | FS 3348 |
|  60 mm | Straight Cutting | 70 cm | FS 3390 | FS 3389 | |
|  90 mm | 3/8 Circle Cutting | 100 cm | | | FS 3397 |

PRE-CUT SUTURES

| SUTURE LENGTH | PACKAGING | SUTURE SIZE IN USP | | |
|---------------|--|--------------------|----------|--------|
| | | 2-0(3) | 1-0(3.5) | 1(4) |
| 5 x 35 cm | Sterilised Monofilament in sterilised foil over wrap. 5 x 35 cm box of 12's pack | FS 703 | FS 704 | FS 705 |
| 2 x 76 cm | Sterilised Monofilament in sterilised foil over wrap. 2 x 76 cm box of 12's pack | FS 903 | FS 904 | FS 905 |