

# DISPOSABLE SKIN STAPLER



*Connecting Tissues*



## Features

- US FDA, CE, ISO and KGMP approved product
- Sterile, ready-to-use, 35 pin surgical stapler
- 316 grade stapler pins to ensure soft but sharp pins
- Unique ratchet system to prevent jamming
- Indicator to check remaining stapler pins
- Ergonomic user friendly design
- Clear view of operating site
- Uniform staple pin formation
- Stapler and staple remover provided at reasonable prices



CE US FDA  
0120 D102093



## FUTURA SURGICARE PVT LTD

No. 29, Balanjaneya Temple Street, Opp. M.S. Ramaiah Hospital,  
KGE Layout, RMV 2nd Stage, Bangalore - 560 094. India

Tel: +91-80-23514126, Fax: +91-80-23514334 E-mail : [futura.surgicare@gmail.com](mailto:futura.surgicare@gmail.com)  
[www.dolphinsutures.com](http://www.dolphinsutures.com)