


Dolphin mesh is available in popular sizes of 2x2, 3x6, 6x6 and 12x12 inches. To facilitate the Surgeons convenience, Dolphin mesh is made available in two types, namely PM and PS types. PM type mesh is thinner, smoother and more flexible. Whereas PS mesh is Strong, thicker and normal flexibility.

Dolphin Meshes are ultrasonically cleaned and packed in internationally acceptable standard pouches, with EO sterile. Further all edges of mesh are laser cut to avoid post surgical inconvenience.

Dolphin Meshes are provided with very High Bursting Strength, Elasticity and tested for Bio-compatibility after implantation. Further the Dolphin Meshes are made available to medical fraternity at reasonable prices.

PROPERTIES	DOLPHIN MESH	
	PM TYPE	PS TYPE
Tensile Strength -cross [in Newtons] -long	220	360
Breaking Elongation % - Cross - Long	120	120
Mullen Bursting Strength [in Newtons]	360	700
Thickness	0.51 mm	0.7 mm
Fabric Weight [gms/m <sup>2</sup> ]	65	120
Fabric Count [dtex]	186	186
Monofilament Diameter	0.15mm	0.15 mm
Pore size	0.85 mm	0.80mm



HK-M

# HERNIA KIT

*Inguinal Hernia Repair*



STERILE EO 

**FS 2407**      **76 cm**  
 Absorbable Surgical Suture U.S.P.  
 Sterilised Surgical Needled Suture.  
**3-0** U.S.P.      **PETCRYL**  
 Undyed  
 Braided Coated polyglycolic Acid (PGA)

**26mm**  
 3/8 Circle Reverse  
 Cutting  


**FS 842**      **90 cm**  
 Non Absorbable Surgical Suture U.S.P.  
 Sterilised Surgical Needled Suture.  
**1-0** U.S.P.      **DURACARE**  
 Monofilament Polypropylene

**30mm**  
 1/2 Circle Round  
 Bodied  


**FS 841**      **90 cm**  
 Non Absorbable Surgical Suture U.S.P.  
 Sterilised Surgical Needled Suture.  
**2-0** U.S.P.      **DURACARE**  
 Monofilament Polypropylene

**30mm**  
 1/2 Circle Round  
 Bodied  


  
 Monofilament Polypropylene (Undyed)  
 Non Absorbable Synthetic Knitted Surgical Mesh

**3.2" x 6"**  
 8cm x 15cm