

NONABSORBABLE SURGICAL SUTURES

DURACARE MONOFILAMENT POLYPROPYLENE [BLUE]

Duracare is a special range of finer monofilament polypropylene sutures made for usage in cardiovascular, ophthalmic, plastic and micro surgeries.

Virtues : Polypropylene sutures are exceptionally smooth and have uniform surface which allows easy passage through tissues. High tensile strength and knot security makes it ideal for cardiovascular usage.



















Material : Synthetic, monofilament polypropylene suture. FDA approved blue colour for easy identification. High quality 300 series alloy needles for exceptional performance in terms of sharpness, strength and ductility.

Tissue Reaction : Minimal tissue reaction. Does not support infection. No tissue rupture.

Applications : Cardiovascular, ophthalmic, plastic and micro surgeries.

Choice : Single and double needle. Available with PTFE pledgets sized 3x3mm or 3x6mm and also without pledgets.

Caution : Avoid damage to the surface of suture

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE IN USP			
				8-0	7-0	6-0	5-0
	3/8 Circle Taper Point		60 cm	CV 2086			
	3/8 Circle Taper Point		60 cm	CV 2088	CV 2078		
	3/8 Circle Taper Point CC		75 cm	CV 8703		CV 8704	
	3/8 Circle Taper Point VB		75 cm	CV 8303	CV 8304	CV 8305	
	3/8 Circle Taper Point		75 cm		CV 8702	CV 8712	
	3/8 Circle Taper Point		60 cm	CV 2080		CV 8804	
	1/2 Circle Taper Point		60 cm				CV 8803
	3/8 Circle Taper Point		75 cm			CV 8802	
	220 micron		75 cm			CV 8597	
	282 micron		75 cm				

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE IN USP				
				6-0	5-0	4-0	3-0	2-0
	3/8 Circle Taper Point		60 cm	CV 8726				
			75 cm	CV 8707				
	3/8 Circle Taper Point CC		75 cm	CV 8706				
	1/2 Circle Taper Point		60 cm			CV 8204		
			90 cm		CV 8710			
	3/8 Circle Taper Point		90 cm		CV 8721			
	3/8 Circle Taper Point VB		75 cm	CV 8307				
	1/2 Circle Taper Cutting		70 cm			CV 899		
	1/2 Circle Taper Cutting		90 cm			CV 8935	CV 8936	CV 8937
	1/2 Circle Taper Point		90 cm		CV 8556	CV 8557		
	1/2 Circle Taper Cutting		90 cm		CV 8556TC			
	1/2 Circle Taper Point		90 cm		CV 8556SA			
	1/2 Circle Taper Point		75 cm		CV 892			
	1/2 Circle Round Bodied		70 cm				CV 822	
	1/2 Circle Taper Cutting		90 cm					CV 877
	1/2 Circle Taper Point		90 cm				CV 8522	
	1/2 Circle Round Bodied		90 cm				CV 825	