











NON-ABSORBABLE SURGICAL SUTURES

PROCARE * Polyester Sutures













- Virtues** : Good tensile strength for long time in tissue. Uniformly coated silicone Provides easy passage through tissues.
- Material** : Precisely braided polyester with silicone coating. Available in green and white colours.
- Tissue Reaction:** Minimal tissue reaction
- Applications** : General closure, Cardiovascular, Plastic surgery, Orthopaedics, Ophthalmic and Skin surgery.
- Choice** : Sterile Needled suture conforming to U.S.P. standards

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE IN USP		
				4-0	3-0	2-0
	1/2 Circle Taper Cutting - Green with PTFE Pledgets of 6mmx3mmx1.5mm		90 cm			PTS 6937
	1/2 circle Taper Cutting - White with PTFE Pledgets of 6mmx3mmx1.5mm		90 cm			PTS 6917
	1/2 Circle Taper Cutting - Green		90 cm	TS 6935	TS 6936	TS 6937
	1/2 Circle Round Bodied - Green		90 cm		FS 623	
	1/2 Circle Round Bodied - Green		100 cm		FS 655	







PROCARE
POLYESTER





NON-ABSORBABLE SURGICAL SUTURES

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE IN USP			
				3-0	2-0	2	5
 25 mm	1/2Circle Taper Cutting White with PTFE Pledgets of 6mmx3mmx1.5mm		90cm		PTS 6987		
 25 mm	1/2circle Taper Cutting Green with PTFE Pledgets of 6mmx3mmx1.5mm		90 cm		PTS 6977		
 25 mm	1/2Circle Taper Cutting Green		90 cm		TS 6977		
 26 mm	1/2Circle Taper Point		100 cm	TS 6552			
 45 mm Heavy	1/2Circle Taper Cutting Heavy		100 cm			TS 4842	
 55mm	1/2Circle T.car Point Heavy		100 cm				TS 4845

SUTURE BOOK

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE IN USP
				2-0
 17 mm	½ Circle Taper Cutting SX 5 Green sutures - 5 White sutures		75 cm	10B52
 26 mm	½ Circle Taper Cutting SX 5 Green sutures - 5 White sutures		75 cm	10B72

SUTURE BOOK WITH PLEDGETS

NEEDLE LENGTH	NEEDLE DESCRIPTION	CUT	SUTURE LENGTH	SUTURE SIZE IN USP
				2-0
 17 mm	1/2Circle Taper Cutting SX Green sutures 5x75 cm - White 5x75cm with PTFE Pledgets 6mmx3mmx1.5mm		75 cm	10B55
 26 mm	1/2Circle Taper Cutting SX Green sutures 5x75 cm - White 5x75cm with PTFE Pledgets 6mmx3mmx1.5mm		75 cm	10B77